




15 Minute *Cryptosporidium* Antibody

NEW - Celllabs *Cryptosporidium* IFA 15 Minute stain

High-definition monoclonal diagnostic test for the rapid detection of *Cryptosporidium* oocysts in human faecal samples as well as water and sewage samples.

Celllabs ... the Australian specialist in waterborne infections are pleased to introduce this high performance diagnostic kit which enables the rapid detection of *Cryptosporidium* oocysts. The monoclonal antibody reacts specifically with the oocyst cell wall to give a clearly defined morphology, making accurate diagnosis possible.



Main features

- Rapid simple assay that detects *Cryptosporidium* oocysts directly
- Positive, reliable results in 15 MINUTES
- Monoclonal antibody specificity
- Identifies environmental sources of infection
- Economical
- Approved by DWI
- **NO NEED FOR REVALIDATION**



Botolph Claydon, Buckingham, MK18 2LR, England

Tel: +44 (0) 1296 714222 ~ Fax: +44 (0) 1296 714806

E-mail: watersciences@tcsgroup.co.uk